

Welcome to ParticleScout



Welcome to the Particle Scout Help.

Installation	Installation Infos
Menu	The ParticleScout Menu
Find Particles	Find and filter particles on (Stitching) Video Images
Particle Manager	Browse and manage all particles
Raman Measurement	Perform Raman measurements on particles
Material Search	Use TrueMatch to perform a material search for each particle spectrum
Report	Create a report and perform statistical analysis

Press the **F1 key** anywhere in the software to open the context help or browse the Help Menu to open the help contents

Installation

Installation

WITec ParticleScout and WITec TrueMatch are both included in the same executable and distributed with the same setup file.

Program Start

ParticleScout is started automatically when

- measuring a Particle Stitching image from the [WITec Control Video Measurement Dialog](#)
- when exporting images from the [Project Manager Context Menu of WITecProject](#).

You can also double-click a .witscout (WITec ParticleScout) file or start the ParticleScout software using the start menu of windows or the desktop shortcut and load a particle project.

Please note that only one instance of WITec TrueMatch / ParticleScout is allowed.


Licensing

A special ParticleScout license is needed to use ParticleScout without any limitations.

Once you have ordered the ParticleScout license, it is included in the WITecProject license and can be used on any number of computers.

Menu

File Menu



Load Particle Project

Loads a previously saved particle project.

Save Particle Project (As)



Saves the current particle project including all particles with their properties, thumbnail images, spectra and database search results.

Exit

Exits the application.

Recent Files

Recent particle scout files are shown at the bottom. Just click it to load the particle project.

ParticleScout Menu**Export Particles****Export Particles to CSV**


Exports all flagged particles with properties into a semicolon separated file.

Export Particle Spectra to ASCII

Exports the measured spectra of all flagged particles into an ASCII file, the format is the [WITec TrueMatch ASCII Import](#) format

Export Selected Particle Spectrum to WITec Project/Control

Exports the measured spectrum of the currently selected particle to the clipboard. Can be pasted as new single spectrum object into WITec Project/Control project manager. Works only if a single particle is selected.

Measurement Options**Offset Correction**

If set, ParticleScout will move the sample positioner with an offset for each particle or if the move to particle button is used.

Hysteresis Compensation

If set, a hysteresis compensation is performed before moving to a particle.

Use Cosmic Ray Removal

If checked, cosmic rays in measured Raman spectra will automatically be removed using a cosmic ray detection algorithm. Only works if at least 2 accumulations are used.

Load Particle Mask from File ...

Loads a bitmap and uses all white pixels as mask pixels for adding new particles. If there are already particles in the current project, you can decide whether to create a new project or add the particles to the current project.

Load Particle Image ...

Loads one or multiple particle image files (e.g. *.bmp or *.png).

If there is already a particle project opened, a message box will open and the user can choose

- to use these images for adding more particles
- to use the image as a main image (if only one image is loaded)
- to create a new project (the first image will be used as the main particle image)

You can load the image with a coordinate information that contains the size, offset and the rotation.

It can be defined by saving a text file (.txt) with the same name as the image file name in the same directory.

Use "." as decimal separator.

Width/Height unit is microns.

X/Y/Z coordinates is the upper left corner of the image, in microns.

Rotation unit is radian.

Pixels must be quadratic.

Example:

Width: 558.24


Height: 558.24

X: 0

Y: 0

Z: 0

Rotation: 0



Save Particle Image ...

Saves the current particle image into a windows bitmap file.

Note that the software allows to save bitmaps with a maximum size of $\approx 700.000.000$ Pixels, which might not be readable in every software.


Particle Image to Clipboard (WIP)

Copies the particle image to the clipboard in WIP format, so it can be pasted into the Project Manager of a running WITecProject or WITecControl instance.

The image will be downsized to 8000x8000 pixels, if larger.

Find Particles

Find Particles



Bright / Dark Particles (Automatically set)

Bright Particles: finds particles that are brighter than the background (e.g. dark field image)

Dark Particles: finds particles that are darker than the background (e.g. bright field image)

Mask by Threshold

Defines which brightness level should be used to detect the particles.

Auto

Calculates an automatic threshold (also used on startup).


For the automatic threshold, "Otsu's method" is used:

https://en.wikipedia.org/wiki/Otsu%27s_method

Mask by Watershed

This algorithm can be used to solve two problems in particle detection.

- Bright and less bright particles are both masked with an individual threshold depending on the particle brightness.
- Particles that are very near to each other can be separated.



To get the best results, please define the parameters step by step in the following order:

Background Threshold

Define a threshold for the background.

When changing this parameter, the preview will temporarily show the background mask.

Detection Threshold

Only areas that are above this threshold (background threshold + detection threshold) will be used to find the particles. When changing this parameter, the preview will temporarily show the detection mask.

Link

If the link check box is checked, the separation and mask thresholds are set equal. This is useful if you are only interested in having individual thresholds for each particle.

Separation Threshold

During the calculation the particles are expanded around local intensity maxima. If the intensity drops below the separation threshold the region becomes a particle by its own. It will not be joined with other particles, which are "falling" below their separation threshold.

Mask Threshold

This parameter defines the mask threshold relative to the local maxima.

Filter Particles

Edge Particles	
<input checked="" type="checkbox"/>	Exclude Edge Particles
Filter Expression	 Quick Filters...
Area > 5  	
Results	
Particles Found:	242
Removed Particles:	74
Particles Used:	168
Define Material	
User Defined Material	PMMA

Exclude Edge Particles

If set, all particles whose contours are not completely within the image are skipped.

Filter Expression / Quick filters

You can filter particles that do not match a custom filter expression.

E.g. you can define that only particles with an area larger than 50 μm should be detected.

See [Filter Expression Editor](#).

User Defined Material

It's possible to export particle images from WITec Project with known material. In this case you can just enter the material name which is then assigned to the particle objects in the resulting particle list.

Actions**Use Results**

This will accept the current mask and proceed with the particle filter.

Preview

The preview image shows the original particle image and the current mask as an overlay.

In the upper right corner you can open options in order to change the brightness / contrast and change the opacity and color of the mask.

Particle Manager

The Particle Manager shows all particles of the current particle project.

It can be used to select particles and use certain particles for Raman measurement, database search or statistical analysis.

Selection for Analysis

Flags

Only flagged particles are used for

- Raman Measurement
- TrueMatch Database Search
- Statistical Analysis
- Export

Selection

The selection (blue rows) can be used to

- flag/unflag particles
- delete particles, particle spectra or particle material
- calculate spectrum properties
- define custom material



Toggle Flag State of Selected Particles

Toggles the flag state of all selected particles.

Select All

Selects all particles in the list.

Select All above Current

Selects all particles in the list that are above the currently selected particle.

Select All below Current

Selects all particles in the list that are below the currently selected particle.

Select via List

You can click on any row to select a single particle. The selected particle will be shown in the Particle Detail View. Use the Shift or Control keys on the keyboard in order to select multiple items.

Toggle Flag in List

You can click the box on the very left side of a particle row entry to flag or unflag this particle. Use the shortcut space to toggle the flags of all selected particles.

Number of Particles

54: Number of flagged particles.

1: Number of selected particles.

55: Total number of particles in the current project.

48: Number of visible particles in the list view.

Particle Filtering

Here you can define a custom filter formula in order to show only particles in the list that match a certain condition.

With the flags buttons you can

- Flag only currently visible particles
- Add currently visible particles to existing flagged items
- Remove currently visible particles from existing flagged items

See [Filter Expression Editor](#).


Sorting

Click on the button in the column header area in order to define which particle property should be shown in the column. The selected property is also used for sorting the items.



If checked, the order of the current list of particles will be reversed.


If particles are sorted by material or any Boolean property, you can quickly jump to the next material / different value using the arrow buttons:

**List View Settings**

Here you can switch between the list view and the thumbnail view. Both views have a custom thumbnail size.

Particle Details

See [ParticleScout Particle Details](#).

Preview Image

The preview image shows the particle image.

All particles, that are visible in the list view are shown as a green overlay.

Red particles are flagged. Blue particles are selected, Pink particles are selected and flagged.

Click to Select


Click on any particle in order to select it in the list. Hold down the left mouse button to change the selection while moving the mouse.



When enabled, click somewhere in the image in order to move the sample positioner to the desired position.

Actions**Find Particles**

Returns to the Find Particles view so you can e.g. change the mask and create a new particle list.

Actions on Flagged Particles**Measurement**

Opens the [Measurement View](#). Here you can acquire a Raman spectrum at each particle and assign it to the particles.

Only available if WITec Control is running.

Material Search with TrueMatch

Opens the [TrueMatch Search View](#). Here you can search all particle spectra in a spectral database and assign a material name to each particle.

Only available if particles with spectra are flagged.

Particle Report

Opens the [Particle Report View](#). Here you create and customize a report of all flagged particles.

Calculate Spectrum Properties

This will calculate an estimation of the amount of Raman or fluorescence signal and a yes-no-info about the over-saturation for all selected particles.

See [Spectrum Properties](#).

Actions on Selected Particles



Define custom material

Here you can enter a material name and assign it to the selected particles.

Delete Spectrum

This will delete the spectrum and its information (material, spectrum properties) from the selected particles. Only available if particles with spectra are selected.


Delete Material

This will delete the material name from the selected particles. Only available if particles with spectra and assigned material are selected.

Delete Particle

Deletes all selected particles.

Particle Details



This view shows the details of a selected or measured particle.

You can select multiple particles in order to see some statistics of the properties of all selected particles: Min, Max, Average, Standard Deviation and Median.



Move to Particle

With this button you can move the sample positioner to the currently selected particle.

Spectrum Viewer

Shows the measured Raman spectrum. See [Spectrum Viewer](#).

Particle Properties

Shows all particle properties.

For a detailed description for each property, see [Particle Properties](#).

Particle Properties

In order to calculate descriptive particle properties from an image the following steps are done:

- The image is converted to a mask
- A list of particles is created (all connected mask pixels belong to one particle)
- For each particle a contour line is calculated
- The pixel positions and the contour line is used to calculate the descriptive particle properties



Definition of Particle Pixels

$$\mathbb{N} = \left\{ \vec{n}_i = \begin{pmatrix} n_x^i \\ n_y^i \end{pmatrix} : 0 \leq i < M_n \right\}$$

\mathbb{N} : Set of all Mask Particles
 M_n : Number of Mask Particles
 \vec{n}_i : Pixel Coordinates

Definition of Particle Contour


$$\mathbb{C} = \left\{ \vec{c}_i = \begin{pmatrix} c_x^i \\ c_y^i \end{pmatrix} : 0 \leq i \leq M_c; \vec{c}_0 = \vec{c}_{M_c} \right\}$$

\mathbb{C} : Set of all Contour Supporting Points
 M_c : Number of Supporting Points
 \vec{c}_i : Pixel Coordinates of Supporting Points

Length and Width

In order to calculate the Length and the Width of the Particle the bounding box with the smallest area is fitted to the contour of the particle.
 The larger side is defined as the Length.

L_{A-Min} : Length of Bounding Box with smallest Area
 W_{A-Min} : Width of Bounding Box with smallest Area
 $L_{A-Min} > W_{A-Min}$



Aspect Ratio

$$\text{Aspect Ratio} = \frac{W_{A-Min}}{L_{A-Min}}$$

The Aspect Ratio range is between 0 to 1.
 A quadratic bounding box has an Aspect Ratio of 1.

Number of Pixels


$$N = M_n$$

N : Number of Particle Pixels

Area


$$A = a \times N$$

A : Particle Area
 a : Pixel Area



Convex Area


$$A_{convex} = A + B$$



Perimeter

$$P = \sum_{i=0}^{M_c-1} |\vec{c}_{i+1} - \vec{c}_i|$$


P : Perimeter



Convex Perimeter


$$P_{convex} = \sum_{i=0}^{M_c-1} |\vec{c}_{i+1} - \vec{c}_i|$$

P_{convex} : Convex Perimeter




Maximum Feret Diameter

The Maximum Feret Diameter is the maximum distance that can be measured if the particle is rotated between a pair of calipers. At least one caliper must touch a segment line of the perimeter.



Minimum Feret Diameter


The Minimum Feret Diameter is the minimum distance that can be measured if the particle is rotated between a pair of calipers. At least one caliper must touch a segment line of the perimeter.



Circular Equivalent Diameter

$$d_{CE} = 2 \sqrt{\frac{A}{\pi}}$$

d_{CE} : Circular Equivalent Diameter




Circularity

$$C = \frac{P_{CE}}{P} = 2 \frac{\sqrt{\pi A}}{P} = \frac{\pi d_{CE}}{P}$$

P_{CE} : Circular Equivalent Perimeter
 C : Circularity

The Circularity range is between 0 and 1.
 If the parameter is 1 the shape of the particle is a circle.



Convexity

$$Convexity = \frac{P_{convex}}{P}$$

The Convexity range is between 0 and 1.
 If the parameter is 1 the shape of the particle is convex.

Solidity


$$Solidity = \frac{A}{A_{convex}}$$

The Solidity range is between 0 and 1.
 If the parameter is 1 the shape of the particle is convex.

Spherical Equivalent Volume

$$V_{SE} = \pi \frac{d_{CE}^3}{6}$$

V_{SE} : Spherical Equivalent Volume



Spectrum Properties

This dialog calculates an estimation of the amount of Raman and Fluorescence signal and a yes-no-info about the over-saturation for all selected particles.

Calculation Options	
Raman and Fluorescence Signal Estimation	
Offset Kind	Range Average
Range for Offset Average [1/cm]	0 0
Noise Threshold Factor	6
Shape Subtraction Size [1/cm]	400
Range for Signal [1/cm]	300 4000
Oversaturation	
Oversaturation Threshold	55000
Range for Oversaturation [1/cm]	100 4000
OK Cancel	

Raman and Fluorescence Signal Estimation

Offset Kind

Before calculating any signal value, a horizontal offset is subtracted from the spectrum:

- Minimum: The minimum value of the spectrum is subtracted from the spectrum
- Range Average: Here you can define a spectral range that is used to calculate an average value (Parameter "Range for Offset Average"). This value is subtracted from the spectrum
- User Defined: A user defined value is subtracted from the spectrum

Range for Offset Average / User Defined Offset

Depending on the selected offset kind, you can define a range for calculating an average value used as offset or directly enter an offset value.

Noise Threshold Factor

This factor defines how much higher the spectral signal must be in comparison to a local noise in order to be detected as Raman signal. The Raman signal value is determined by all spectral pixels that are higher than $\langle \text{Noise Threshold Factor} \rangle * \langle \text{Local Noise} \rangle$.

Shape Subtraction Size

Defines the shape size of the [shape background subtraction](#). The fluorescence signal value is determined by the subtracted background.

Range for Signal

Defines the spectral range that is used for the signal estimation.

Oversaturation

Oversaturation Threshold

If any spectrum pixel value within the defined range is higher than the threshold, the spectrum is marked as over-saturated.

Range for Oversaturation

Defines the spectral range that is used for the over-saturation calculation. The Rayleigh Peak might be saturated, so you can e.g. skip the Rayleigh area.

Raman Measurement

The Raman Measurement view allows to acquire a Raman spectrum for each particle.

WITec Control Configuration Setup

Before starting the Raman Measurement, you can choose one of the Raman configurations in WITec Control (e.g. "Raman CCD 1").

Its possible to measure a set of particles, then switch the WITec Control configuration to "Raman CCD 2" or change to another excitation laser, then measure other particles.

Measurement Sequence

For each particle, the following tasks are done by the software:

- The Sample Positioner is moved to the particle using the [Positioning Settings](#)
- Depending on the Z-Axis Behavior, a spectral Auto Focus is performed
- The spectrum is measured

Single Spectrum

Single Spectrum	
Measurement Mode	Optimize Fast ▾
Accumulations	100
Integration Time [s]	0.1000
Low Signal Limit	50
SNR Limit	40
Edit Mask	

Measurement Mode

- Normal: Uses the defined number of accumulations
- Optimize: Only accumulates spectra that improve signal to noise ratio of the Raman signal. Spectra with high fluorescence signal are rejected.
- Optimize Fast: In addition to the above option this mode will stop before the defined number of accumulations is reached:
 - if the signal to noise limit of the optimal accumulated spectrum is reached.
 - if the Raman signal of a single spectrum is too low
 - if the signal to noise limit of the optimal accumulated spectrum is not improving.

Accumulations

The number of accumulations for the spectrum measurement.

If more than 1 accumulation is used, a cosmic ray removal can be performed.

Integration Time

The integration time for each spectrum accumulation.

Low Signal Limit

Only used for Measurement Mode "Optimize Fast".

If the signal of a single spectrum is lower than this value the measurement will stop in order to save time.

SNR Limit

Only used for Measurement Mode "Optimize Fast".

If the optimal accumulated spectrum has a higher signal to noise ratio the measurement will stop in order to save time.

Edit Mask

Only used for Measurement Modes "Optimize" and "Optimize Fast".

Here you can define a spectral mask to define which parts of the spectrum should be used to calculate the signal.

Z-Axis Behavior

Z-Axis Behavior	
<input type="radio"/> No Z Movement	
<input checked="" type="radio"/> Spectral Autofocus	
<input type="radio"/> Fix Z Position	0.00

- No Z Movement: user can move the Z-Axis to a desired location
- Spectral Autofocus: a spectral auto focus is performed at each particle
- Fix Z Position: lets you define an absolute Z-Position (Software controlled, limited Z-axis space)

Spectral Auto Focus

Before using the spectral Auto Focus adjust the following parameter:

Spectral Auto Focus	
Z-Axis Range [µm]	<input type="text" value="-10.0"/> <input type="text" value="90.0"/>
Min. Integration Time [s]	<input type="text" value="0.0500"/>
Step Size Multiplier	<input type="text" value="1.0"/>
<input type="button" value="Edit Mask"/>	
<input type="button" value="Execute Spectral Auto Focus"/>	

Z-Axis Range / Min. Integration Time / Step Size Multiplier

See [Spectral Auto Focus](#) documentation of WITec Control.

Edit Mask

Lets you edit a spectral mask to define which parts of the spectrum should be used for performing the auto focus.

Execute Spectral Auto Focus

Executes the spectral auto focus at the current position.

This way you can test if the spectral auto-focus works as expected.

Measurement Order and Additional Options

Measurement Order
<input type="button" value="Shortest Path"/>

Additional Options
<input type="button" value="Positioning and Preprocessing"/>

Measurement Order

Here you can define in which order the particles are measured:

- Shortest Path: calculates an estimation of the shortest path in order to save time while traveling to each particle.
- User Defined: uses the sort order of the list view in the Particle Manager, e.g. "Area".

Positioning and Preprocessing

See [Measurement Options](#).

Actions

**Use Results**

This will assign the measured spectra to the particles.

**Start Measurement**

Starts the measurement.

**Pause / Stop Measurement**

Stops the measurement.

You can continue at the last measured particle or measure all particles again.

Progress

Measured Particles 18 of 55	Measurement Running. Remaining Time: 00:02:08
---------------------------------------	--

Shows the progress and number of measured particles.

Preview

Preview
Selected Preview <input type="text" value="13"/> <input type="range"/>

Here you can select which measured spectrum and corresponding particle should be shown.

If the last measured particle is selected, the preview will automatically stay at the latest measurement.

Particle Details

See [Particle Details](#).

Material Search

The spectra of all flagged particles are used in TrueMatch in order to search for chemical components.

The names of the best component search results are assigned to the particles as material name and HQI (hit quality).

This information can be used:

- to sort, filter and select particles in the Particle Manager
- to decide if certain particles spectra should be measured again, e.g. with different measurement parameters
- to create a report with statistics

See [TrueMatch Search Overview](#).

Actions



Use Results

This will assign the best search results as material name to the particles.

Result Selection

Please make your Selection

Selected Result Component 1

Minimum Weight [%] of Selected Component

If a multi-component search was performed, you can define which of the sub-results should be used (Component 1, 2, 3). You can also define a minimum weight percentage for the selected component in order to avoid using results with a very low importance.

Report

The report view allows you to define a report layout with the following elements:

- Particle Image (Map) with Legend
- Typical Particles with Thumbnails and Spectra
- Bar Chart with Legend
- Pie Chart with Legend
- Table with Categories

The resulting report view can then be exported as bitmap into the windows clipboard.

Report Settings

Settings

Configuration 1 

"Configuration 1" ComboBox

Here you can select between your saved configurations.

Options Button


Opens the report configuration window in order to define the report elements and data preparation.

See [Report Configuration](#).

Bitmap to Clipboard

Copies the whole report as a bitmap to the clipboard.

Example:



Report Configuration

The ParticleScout Report is configured using the XML Format ([Extensible Markup Language](#)) with the following structure:

```

WITec
Style
Data
  Colors
  Groups
  Categories
  Visual Elements (Map, Thumbnail, Bars, Pie, Table)
Data
[...]
```

Press on the examples button in the software in order to create example xml code automatically.

The script must begin with
<WITec>

and end with
</WITec>

Style Tag

You can define exactly one Style Tag.

The Style Tag can be used to define some display properties:

```

<Style ScaleFactor='2'
  ShowCaptions='true'
  CaptionBackground='LightGray'
  CaptionFontSize='16'
  CaptionPlacement='Top'
  CaptionBorderThickness='1'
  CaptionBorderColor='Gray'
  ElementBackground='Yellow'
  CanvasBackground='#11000000'>
```

ScaleFactor: Defines a scale factor to scale the size of all report drawings (0.3 to 3, default: 1).

ShowCaptions: Show or hide all captions (default: true).

CaptionBackground: Defines the color of the caption background (default: #11000000).

CaptionFontSize: Defines the caption font size (default: 14).

CaptionPlacement: Can be top or bottom (default: Top).

CaptionBorderThickness: Defines the border size around the caption (default: 0).

CaptionBorderColor: Defines the color of the caption border (default: Transparent).

ElementBackground: Defines the color each elements background (default: Transparent).

CanvasBackground: Defines the color of the canvas (default: White).

Data Tag

You can define any number of data tags. At least one data tag must exist.

Within each data tag you can define colors, sorting, grouping and categorization of particle properties.

Any number of visible elements can be specified in order to define how the sorted data should be presented in the report.

```

<Data GroupCategoryName='Material'
  SortCategory='Total Area'>
  <Colors> [...]
```

```

<DynamicGroup> [...]
<NamedGroup> [...]
  <Category> [...]
  <Map> [...]
  <Bars> [...]
  <Pie> [...]
  <Thumbnail> [...]
  <Table> [...]
</Data>

```

GroupCategoryName: Is only used for tables. Defines the header of the group column.

SortCategory: Defines the category name of the category that should be used for sorting the data.

Colors Tag

Colors are used for maps and bar charts.

The number of defined colors will define the number of categories visible in maps and bar charts.

You can define a color by using hexadecimal (A)RGB Values or by using the string representation (please refer the System.Windows.Media.Colors class in Microsoft documentation):

```

<Colors DefaultColor='LightGray'>
  <Color Value='#4572A7' />
  <Color Value='#4BACC6' />
  <Color Value='#92D050' />
  <Color Value='#F79646' />
</Colors>

```

DefaultColor: Is only used for Table elements. Defines the group name / caption of the most left column.

DynamicGroup Tag

This will automatically create a dynamic number of groups of particles.

Each group is defined by a unique string/name which is evaluated from each particle using a string expression.

```
<DynamicGroup NameExpression='Material' Condition='NOT Material.Equals("PET")' />
```

NameExpression: A string expression.

This expression will be automatically evaluated for each particle.

You can use any particle property here, for example Material.

All particles returning the same string will be in the same group (e.g. all particles with the Material "PMMA").

Condition: A boolean expression (optional). Only particles matching this expression will be used for creating groups.

NamedGroup Tag

This will create a group with a user defined name and condition.

Name: The name of the group. Used as table row group name.

Condition: A boolean expression. Only particles matching the expression will be part of the group.

Color: Color of the group. Used in map and map legend.

Position: Can be "Top", "Bottom" or "Default". Makes it possible to put a group at the top or bottom of all sorted groups.

Category Tag

This will create a category with a user defined name and condition. Categories can be used as columns in table presentation or in bar charts.

Name: The name of the category. Used as table column header.

Condition: A boolean expression. Only particles matching the expression will be part of the category.

Count: Allows to define how the "sum value" of the category is calculated.

A value of "1" is used, if no "Count Expression" is defined.

Must be a method call of one of the following methods:

Min, Max, Sum, Average, Median, StandardDeviation.

Usage Example: "Sum(Area)". The string in brackets is a numeric expression that can use all of the particle properties.

NormMode: Normalizes the category results for each group by dividing by the count result of all particles (None, Relative1, Relative100).

Color: Color of the category. Only used as table column background color.

StringFormat: The string representation expression for table cell content. Please refer to the Microsoft .NET string format documentation.

Example: If the number is 0.2345, then StringFormat "F0" = 0, "F1" = 0.2, "F2" = 0.23, "P2" = 23.45 %

Visual Elements

The following properties can be used in all visual elements:

Caption: Used for a surrounding group box header.

SideBySide: Set to "true" to place the element horizontally right next to the previous visual element.

ShowEmptyGroups: Set to "true" to show groups that have no matching particles.

ShowEmptyCategories: Set to "true" to show categories that have no matching particles.

InvertData: Set to "true" to display groups as categories and vice versa.

Map Tag

This will present the groups in a map using the defined colors.

```

<Map Background='Transparent'
  ImageOpacity='0.0'
  MaskOpacity='1.0'

```

```

Width='400'
LegendType='NameAndValue'
ShowScaleBar='true'
ScaleBarForeground='Blue'
ScaleBarBackground='White'
ScaleBarVerticalAlignment='Bottom'
ScaleBarHorizontalAlignment='Left' />

```

Background: The background color (used if ImageOpacity!=1, default: Transparent).

ImageOpacity: Defines the opacity of the particle image (0 to 1, default: 0).

MaskOpacity: Defines the opacity of the particle mask (0 to 1, default: 1).

Width: The width of the map in pixels (2 to 4000, default: 400).

The height is automatically adjusted using the correct ratio.

The pixel size scales with the global report scale factor.

LegendType: Visibility of the legend (None, Name, NameAndValue, default: NameAndValue).

ValueCategory: The name of a Category that defines the values for the legend (default: not defined).

ValueStringFormat: String format for the values (default: not defined).

ShowScaleBar: Set to "true" to show a scale bar on the map (default: true).

ScaleBarForeground: The scale bar color (default: Blue).

ScaleBarBackground: The scale bar background color (default: White).

ScaleBarVerticalAlignment: The scale bar vertical position (Bottom, Top, Center, default: Bottom).

ScaleBarHorizontalAlignment: The scale bar horizontal position (Left, Center, Right, default: Left).

Bars Tag

Shows a bar chart with an X/Y Axis and legend.

```

<Bars Caption="Categories"
Width='500'
Height='300'
BarsSideBySide='false'
Normalize='true'
BarWidth='0.7'
XAxisLabelRotation='45'
LegendType='NameAndValue'
InvertData='true' />

```

Width: The width of the chart, in DPI depending units (4 to 4000, default: 500).

Height: The height of the chart, in in DPI depending units (4 to 4000, default: 300).

BarsSideBySide: Set to "true" to show each group as a separate bar next to each other (default: false)

Normalize: Set to "true" to normalize categories to 100% (default: false).

BarWidth: The width of each category, as factor (0.1 to 1, default 0.7).

XAxisLabelRotation: The X Axis label angle in degrees (-180 to 180, default: 45)

LegendType: Visibility of the legend (None, Name, NameAndValue)

ValueCategory: The name of a Category that defines the values (default: not defined).

ValueStringFormat: String format for value labels (default: not defined).

Pie Tag

Shows a pie chart with legend.

```

<Pie Caption="Some Pie"
Width='300'
Height='300'
OutsideLabelStyle='NameAndValue'
InsideLabelStyle='None'
LegendType='Name' />

```

Width: The width of the chart, in DPI depending units (4 to 4000, default: 300).

Height: The height of the chart, in in DPI depending units (4 to 4000, default: 300).

OutsideLabelStyle: Visibility of the label outside from each pie segment (None, Name, NameAndValue, Value, default: None).

InsideLabelStyle: Visibility of the label inside each pie segment (None, Name, NameAndValue, Value, default: None).

BorderColor: Color of the border of each pie segment (default: White)

BorderThickness: Thickness of the border of each pie segment, in DPI depending units (default: 1.0)

LegendType: Visibility of the legend (None, Name, NameAndValue, default: NameAndValue)

ValueCategory: The name of a Category that defines the values (default: not defined).

ValueStringFormat: String format for value labels (default: not defined).

Thumbnail Tag

This will present a list of particle groups with name, particle thumbnails and a representing spectrum (if measured).

```

<Thumbnail NumberOfEntries='3'
NumberOfThumbnails='3'
ThumbnailWidth='60'
SpectrumWidth='300'
SpectrumHeight='60'
ShowXAxis='false'
ShowYAxis='false'
XAxisTitle='1/cm'
SortThumbnailsBy='Area'
SortSpectraBy='HQL' />

```

NumberOfEntries: The maximum number of list entries (maximum is number of groups, default: 3).

NumberOfThumbnails: The number of thumbnails (0 to 20, default: 2).

ThumbnailWidth: The width of each thumbnail, in DPI depending units (0 to 4000, default: 60).

SpectrumWidth: The width of the spectrum, in in DPI depending units (0 to 4000, default 400).

SpectrumHeight: The height of the spectrum, in in DPI depending units (0 to 4000, default: 60).

ShowXAxis: Set to "true" to show the X Axis.

ShowYAxis: Set to "true" to show the Y Axis.

XAxisTitle: The title/caption of the X Axis.

SortThumbnailsBy: An expression returning a particle property used for sorting the thumbnails (e.g. "Area" -> the particles with the biggest areas are shown)

SortSpectraBy: An expression returning a particle property used for sorting the spectra (e.g. "Area" -> the spectrum of the particle with the biggest area is shown)

Table Tag

Shows a table with groups as rows and categories as columns.

```
<Table ColumnWidth='50'
  MasterColumnWidth='0'
  HeadingColorOpacity='0.5'
  CellBackground='LightBlue'
  CellBorderThickness='1'
  CellBorderColor='Transparent'
  CellPadding='2'/>
```

ColumnWidth: Defines the width of each category column, in DPI depending units. If 0, the width is minimal for each column (default: 0).

MasterColumnWidth: Defines the width of the most left column, in DPI depending units. If 0, the width is minimal (default: 0).

HeadingColorOpacity: Defines the opacity of the heading colors or columns and row headings (0 to 1, default: 0.5).

CellBackground: Defines the background color for each cell (default: #19000000).

CellBorderThickness: Defines the thickness of each cells border. See [thickness format](#). (default: 1).

CellBorderColor: Defines the color of the cell border. Only visible if the thickness is not 0. (default: Transparent).

CellPadding: Defines the space around cell content. See [thickness format](#). (default: 2).

Thickness Format

Use one number to define the same thickness for left/right/top/bottom: '5'

Use two numbers to define left/right and top/bottom: '5 2'

Use four numbers to define left/top/right/bottom: '2 2 1 0'

Possible Color Strings

	AliceBlue	#FFF0F8FF		DarkTurquoise	#FF00CED1		LightSeaGreen	#FF20B2AA		PapayaWhip	#FFF9EFD5
	AntiqueWhite	#FFFAEBD7		DarkViolet	#FF9400D3		LightSkyBlue	#FF87CEFA		PeachPuff	#FFF9DAB9
	Aqua	#FF00FFFF		DeepPink	#FFFF1493		LightSlateGray	#FF778899		Peru	#FFCD853F
	Aquamarine	#FF7FFFD4		DeepSkyBlue	#FF00BFFF		LightSteelBlue	#FFB0C4DE		Pink	#FFC0C0CB
	Azure	#FFF0FFFF		DimGray	#FF696969		LightYellow	#FFFFFFE0		Plum	#FFDDA0DD
	Beige	#FFF5F5DC		DodgerBlue	#FF1E90FF		Lime	#FF00FF00		PowderBlue	#FFB0E0E6
	Bisque	#FFF7E4C4		Firebrick	#FFB22222		LimeGreen	#FF32CD32		Purple	#FF800080
	Black	#FF000000		FloralWhite	#FFFFFFAF0		Linen	#FFFAF0E6		Red	#FF000000
	BlanchedAlmond	#FFF9EBCD		ForestGreen	#FF228B22		Magenta	#FFFF00FF		RosyBrown	#FFB0C8F8
	Blue	#FF0000FF		Fuchsia	#FFFF00FF		Maroon	#FF800000		RoyalBlue	#FF4169E1
	BlueViolet	#FF8A2BE2		Gainsboro	#FFDCDCDC		MediumAquamarine	#FF66CDAA		SaddleBrown	#FF8B4513
	Brown	#FFA52A2A		GhostWhite	#FFF8F8F8		MediumBlue	#FF0000CD		Salmon	#FFFA8072
	BurlyWood	#FFDEB887		Gold	#FFFD7000		MediumOrchid	#FFBA55D3		SandyBrown	#FFD4A460
	CadetBlue	#FF5F9EA0		Goldenrod	#FFDAA520		MediumPurple	#FF9370DB		SeaGreen	#FF20B2B7
	Chartreuse	#FF7FFF00		Gray	#FF808080		MediumSeaGreen	#FF3CB371		SeaShell	#FFF5F5EE
	Chocolate	#FFD2691E		Green	#FF008000		MediumSlateBlue	#FF7B68EE		Sienna	#FFA0522D
	Coral	#FFF77F50		GreenYellow	#FFADDF2F		MediumSpringGreen	#FF00FA9A		Silver	#FFC0C0C0
	CornflowerBlue	#FF6495ED		Honeydew	#FFF0FFFF		MediumTurquoise	#FF48D1CC		SkyBlue	#FF87CEEB
	Cornsilk	#FFFFFF8DC		HotPink	#FFF69B4		MediumVioletRed	#FFC71585		SlateBlue	#FF6A5ACD
	Crimson	#FFDC143C		IndianRed	#FFCD5C5C		MidnightBlue	#FF191970		SlateGray	#FF708090
	Cyan	#FF00FFFF		Indigo	#FF4B0082		MintCream	#FFF5FFFA		Snow	#FFF5F5FA
	DarkBlue	#FF00008B		Ivory	#FFFFFFF0		MistyRose	#FFF5E4E1		SpringGreen	#FF00FF7F
	DarkCyan	#FF008B8B		Khaki	#FFF0E68C		Moccasin	#FFF5E4B5		SteelBlue	#FF4682B4
	DarkGoldenrod	#FF8B860B		Lavender	#FFE6E6FA		NavajoWhite	#FFF5DEAD		Tan	#FFD2B48C
	DarkGray	#FFA9A9A9		LavenderBlush	#FFF5F0F5		Navy	#FF000080		Teal	#FF008080
	DarkGreen	#FF006400		LawnGreen	#FF7CFC00		OldLace	#FFF5F5E6		Thistle	#FFD8BFD8
	DarkKhaki	#FFBDB76B		LemonChiffon	#FFFFACD0		Olive	#FF808000		Tomato	#FFF56347
	DarkMagenta	#FF8B008B		LightBlue	#FFADD8E6		OliveDrab	#FF6B8E23		Transparent	#00FFFFFF
	DarkOliveGreen	#FF556B2F		LightCoral	#FFF08080		Orange	#FFFA5000		Turquoise	#FF40E0D0
	DarkOrange	#FF8F0000		LightCyan	#FFE0FFFF		OrangeRed	#FF450000		Violet	#FFEE82EE
	DarkOrchid	#FF9932CC		LightGoldenrodYellow	#FFFAFAD2		Orchid	#FFDA70D6		Wheat	#FFF5DEB3
	DarkRed	#FF8B0000		LightGray	#FFD3D3D3		PaleGoldenrod	#FFEE8AA8		White	#FFFFFFF5
	DarkSalmon	#FFE9967A		LightGreen	#FF90EE90		PaleGreen	#FF98FB98		WhiteSmoke	#FFF5F5F5
	DarkSeaGreen	#FF8FBC8F		LightPink	#FFF5B6C1		PaleTurquoise	#FFAFEEEE		Yellow	#FFFFF000
	DarkSlateBlue	#FF483D8B		LightSalmon	#FFFA07A		PaleVioletRed	#FFDB7093		YellowGreen	#FF9ACD32
	DarkSlateGray	#FF2F4F4F									

Filter Expression Editor

The Filter Expression Editor allows to filter particles before [creating the particle list](#) or to hide or select certain particles in the [Particle Manager](#).


Formula Editor

```
(Area > 50 and NumPixels > 20)
or
Circularity > 0.5
or
Area < 3
```

Enter a boolean formula that defines a filter condition.
Click on the yellow question mark in order to see all possible variables and some examples.


To remove the formula, press X.

If there is some error, you can hover or click the red exclamation mark in order to see what's wrong:



Quick Filters


You can define quick filters in order to recall a custom formula with a single mouse-click.



Click on a desired Quick Filter to use the saved formula in the filter expression editor.

In Particle Manager: click on the desired Quick Filter check mark to add/remove particle flags.

Configure Quick Filters



Here you can add, remove or edit Quick Filters.

Just add a new filter and set a name and formula.

Example Formulas

Area > 5	Larger than condition
Area > 5 and Area < 10	Combination with AND
(Area > 10 and Area < 20) or (Area > 50 and Area < 70)	Combination with AND and OR using brackets
Material = "Quartz"	String comparison
Material.Contains("Quar")	String method call "Contains"
FeretMin > 3 * FeretMax	Larger than condition with multiplication
RandomValue > 0.5	Random selection of half of all elements (RandomValue = Random number between 0 and 1)
IsOversaturated	Show only particles that have over-saturated spectrum
!IsOversaturated	NOT operator "!": show only particles that have NO over-saturated spectrum
Math.Sqrt(Length * Width) > 5	Calculates the square root of Length * Width and compares larger than 5.

Please refer to the .NET Framework 4.7.1 documentation for other Math or String methods.

